

Bedside cabinet with overbed table – bilateral functionality

Crafted to provide a welcoming home environment in today's care facilities.



B15



Features and benefits

- Drawers that glide effortlessly into place
- Long-lasting, durable construction, using materials that are both easy to clean and maintain
- Aptly sized to suit the needs of today's care environments



Optional: Square top plate



Optional: Top railing



Optional: Plastic drawer insert



Optional: Towel holder



Optional: 2 x single bottle holder


Optional: Contrasting colours front
(drawers and doors)

Optional: Contrasting colours front
and overbed table

Technical Specifications model SER - CH775A

Height	800 mm
Width	600 mm
Depth	460 mm
Weight	34.5 kg
Overbed table	
Height adjustment	610 to 950 mm
Surface area	547 x 349 mm
Tilting range	max. 60°
Material	
Body / Top plate	Melamine, Laminate
Edges	Solid wood
Overbed table	High Pressure Laminate (HPL)
Castors	
	5 x 65 mm with 2 brakes

Top plate	
standard	round
optional	square
Handles/Colours	
see handles and colours sheet	
Additional options	
<ul style="list-style-type: none"> • top railing • plastic drawer insert • towel holder • x 2 single bottle holder • contrasting colours for <ul style="list-style-type: none"> - body / top plate and /or overbed table - front (drawers and doors) and /or overbed table - overbed table 	

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries.
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: SEN150328-01, 25 November 2015