

# Bedside cabinet with overbed table - BL260 bilateral functionality

Offering practical design for today's patient environment



## Features and benefits

- Ideal for all budgets, with a wide variety of colour options
- Long-lasting, durable construction, using materials that are easy to both clean and maintain
- Large storage capacity to suit the needs of today's patient environments

## Technical Specifications **model BL260**

### Height

with 50 mm castors	963 mm
with 75 mm castors	985 mm

<b>Width</b>	559 mm
--------------	--------

<b>Depth</b>	464 mm
--------------	--------

<b>Weight</b>	~39 kg
---------------	--------

### Overbed table

height adjustment	711 to 1,111 mm
surface area	340 x 530 mm
tilting range	360°

### Material

top plate (removable)	plastic
body	laminate
edges	aluminium
base	polyurethane
handles	metal

### Castors

standard	4 x 50 mm double-band with 2 brakes
optional	4 x 75 mm double-band with 2 brakes

### Colours

see colours sheet

### Additional options

2 accessory bar holders

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health**. Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers**.

Not all products/options are available in all countries.  
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

**[www.hill-rom.com](http://www.hill-rom.com)**

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

#### Furniture manufacturer:

Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN150366-02, 07 April 2016