

Overbed table – RT-985

Sturdy functionality with ergonomic design

Unverwechselbar
JÜLKER™



RT-985 (tilting function)

Features and benefits

- Efficient U-shaped undercarriage means the table is compatible with all care environments, delivering maximum working stability
- Ergonomic design features a reduced-width front wheelbase, allowing it to be positioned in all areas of the patient/resident environment
- Locking dual-directional tilted position offers patient/resident an ideal reading position
- The fully-rotational table top and aluminium undercarriage facilitate easy cleaning

Overbed table – RT-985

Technical Specifications **model RT-985**

Height adjustment (without rails)	740 – 1,165 mm
--	----------------

Width	975 mm
--------------	--------

Depth	450 mm
--------------	--------

Surface area

standard	tilting portion: 675 x 420 mm fixed portion: 210 x 420 mm
optional	fixed table: 885 x 420 mm

Weight	~20 kg
---------------	--------

Maximum load	15 kg
---------------------	-------

Material

surface area	High Pressure Laminate (HPL)
undercarriage	Aluminium

Castors

standard	4 x 65 mm double-band with 2 brakes
----------	---

optional	4 x 50 mm or 4 x 75 mm double-band with 2 brakes
----------	---

Tilting range

360° rotational

Colours

see colours sheet

Hill-Rom is a leading global medical technology company with more than 10,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries.
For further information about our products and services, please
contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN150333-01, 07 October 2015

Enhancing outcomes for
patients and their caregivers:

Hill-Rom