



# Bariatric Accessories

Comprehensive Bariatric Care

Comprehensive solutions to help safely streamline routine tasks and provide a comfortable experience for both patients and caregivers.

Hospitals are adjusting to accommodate an increasing population of bariatric patients. To help minimize the risks associated with manual handling of obese individuals, staff must be aware of what equipment is available, how to access it quickly and trained to use with bariatric patients. Unfortunately, many healthcare facilities lack equipment and protocols for lifting and moving bariatric patients.

Key Features at a Glance:

- **Bariatric Expandable Walker** — Adjustable width and height to accommodate various patient sizes.
- **Bariatric Wheelchair** — Swing-away elevating leg rests allow for easier patient ingress and egress.
- **Extra Large Commode** — Easy-to-release drop arm allows you to make safe lateral transfers.
- **Shower Transfer Bench** — Durable, molded plastic bench and backrest provides adjustable height of individual zones.

Hill-Rom can help **transform bariatric care.**

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**

## Bariatric Accessories

### Comprehensive Bariatric Care



#### Technical Specifications

##### Wheelchair\*

Seat width	30" (76.2 cm)
Seat depth	20" (50.8 cm)
Maximum weight capacity	700 lbs (317 kgs)

##### Expandable Walker\*

Height	Adjusts 29"-36" (73.7-91.4 cm)
Width inside handle grip	19"-25" (48.3-63.5 cm)
Depth opened at base	25" (63.5 cm)
Maximum weight capacity	1,000 lbs (453 kgs)

##### Extra Large Commode\*

Seat height	18.5" (47 cm)
Seat width	36" (91.4 cm)
Seat depth	21" (53.3 cm)
Maximum weight capacity	1,000 lbs (453 kgs)

##### Shower Transfer Bench\*

Seat height	18.25"-23.25" (46.35-59.05 cm)
Seat width	27" (68.6 cm)
Seat depth	16" (40.6 cm)
Maximum weight capacity	600 lbs (272 kgs)

\* Produced by Drive Medical Design and Manufacturing

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
CTG128 rev 8 29-OCT-2015 ENG - US

For further information about this product or a service,  
please contact your local Hill-Rom representative or  
visit our webpage:

**USA 800-445-3730**  
**Canada 800-267-2337**

**[www.hill-rom.com](http://www.hill-rom.com)**

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**