

3080

Long Term Care bed for maximum versatility and stability

Unverwechselbar
VÖLKER™



Features and benefits

Telescopic split siderails provide **maximum protection** during articulation of the bed

Easy-to-grip split siderails provide **secure support** to enable safe and independent egress, which encourages resident independence

When elevated, the **telescopic split siderails outline the perimeter** of the bed to restore residents' feeling of security

Residents can **independently adjust** height, lying surface and sitting position using the intuitive, ergonomically designed pendant control

The innovative **Völker Micro-stimulation System MiS®** can lead to enhanced residents' comfort and improved sleep

The stable, four-column undercarriage ensures a **safe working position** at any height

Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Völker 3080



Standard: Pendant control



Optional: Caregiver pendant control

Technical Specifications

Safe Working Load	230 kg
Weight	96 kg
Height adjustment	
standard	40-80 cm
optional	35-70 cm
Protection class	IPX4
Bed angles	
Head section	0-70°
Knee section	0-45°
Reverse Trendelenburg	+14°
Sleep deck width	
standard	90 cm
optional	100 cm
Sleep deck length	
standard	200 cm
optional	210 or 220 cm
Sleep deck surface (removeable)	
standard	Micro-stimulation System MiS®
optional	High Pressure Laminate (HPL)

Split siderail height (dependant on sleep deck surface)

standard	34-35.5 cm
optional	37-38.5 cm
optional	43.5-45 cm

Castors	35 mm single 100 kg/castor
---------	-------------------------------

Lifting pole holder	inside double holder
---------------------	----------------------

Pendant controls	
standard	standard
options see pendant controls sheet	

Caregiver pendant controls	
options see pendant controls sheet	

Designs/Colours	
see Völker Long Term Care beds brochure	

Additional options

Linen holder
Manual bilateral CPR
Protective wall bumpers

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries.
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN139303-02, 09 February 2016

Enhancing outcomes for
patients and their caregivers:

Hill-Rom