

5384 Kepler™ low-height bed - Style Unverwechselbar DÜLKER™

The evolution in Long Term Care



Design S

Features and benefits

The low minimum height of 20 cm* helps **reduce the severity of injury** should a fall occur, and increases the range of positioning for caregivers

Patented Bed Height Memorisation allows residents and caregivers to easily select a preferred height and position

The ability to raise the bed to 80 cm* offers an **optimal ergonomic working position** for caregivers

Residents can **independently adjust** height, lying and sitting position using the intuitive pendant control

Automatic egress positioning contributes to safer and easier bed exit

The **aesthetic full-length siderails** are integrated into the bed when lowered

"Safe low" positioning exclusively for caregiver use assists in the prevention of unforeseen entrapment

Automatic level and horizontal positioning enhances ease-of-use for the caregiver

** dependent on configuration chosen*

Enhancing outcomes for
patients and their caregivers:

Hill-Rom

5384 Kepler low-height bed - Style



Optional: linen holder



Optional: 4 manoeuvrable castors



Optional: EasyChair™ positioning



Optional: MiS® sleep deck surface

Technical Specifications (standard configuration)

Safe Working Load	230 kg
Weight	130 to 159 kg*
Height adjustment	20 to 80 cm*
Height memorisation range	35 to 55 cm* or 40 to 60 cm*
Protection class	IPX4
Bed angles	
Head section	~ 0 to 70°
Knee section	~ 0 to 45°
Reverse Trendelenburg	+14°
Sleep deck width	90 cm
Sleep deck length	200 cm
Sleep deck surface (removable)	
standard	polymer
optional	Micro-stimulation System MiS®
Full-length siderail height	34 to 35.5 cm*
Lifting pole holder	inside double holder (adapter required for IV pole)

* dependent on configuration chosen

Castors

standard	
head end (fixed)	2 x 80 mm single-band
foot end (manoeuvrable)	4 x 100 mm double-band
optional	
head end (manoeuvrable)	2 x 100 mm double-band
foot end (manoeuvrable)	2 x 100 mm double-band

Pendant controls, designs and colours

see pendant controls, designs and colours sheet

Manual bilateral CPR

Additional options

Linen holder
Protective wall bumpers
Backup battery
Secondary pendant control
Bilateral accessory bars including restraint mounts
EasyChair positioning**
Night light**
Trendelenburg**

** available only as a bundle

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries.
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN139308-02, 15 February 2016

Enhancing outcomes for
patients and their caregivers:

Hill-Rom