

Vis-a-Vis

Advanced mobility for participation in daily activities

Unverwechselbar
VÖLKER™



Features and benefits

The sitting position is effortlessly achieved enabling **face-to-face communication** between residents and caregivers/visitors

Telescopic split siderails provide **maximum protection** during mobilisation

Easy-to-grip split siderails provide **secure support** to enable safe and independent egress, which encourages resident independence

The innovative **Völker Micro-stimulation System MiS®** can lead to enhanced resident comfort and improved sleep

The advanced design enables individual caregivers to mobilise residents to a sitting position safely and easily, **reducing the risk of caregiver injuries** resulting from resident repositioning

The unique healthcare design from Völker ensures that beds give the resident's room a **home-like ambience** to promote wellbeing

Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Völker Vis-a-Vis



Autocontour Comfort pendant control with siderail clip (horizontal)



Caregiver pendant control with total lock-outs



Optional: Sleep deck surface (HPL)

Technical Specifications

Safe Working Load	250 kg
Weight	169 kg
Height adjustment	40-70 cm
Protection class	IPX4
Bed angles	
Head section	0-70°
Knee section	0-45°
Reverse Trendelenburg/Trendelenburg	+14°/-14°
Sleep deck width	
standard	90 cm
optional	100 cm
Sleep deck length	210 cm
Sleep deck surface (removeable)	
standard	Micro-stimulation System MiS®
optional	High Pressure Laminate (HPL)
Split siderail height (dependant on sleep deck surface)	
standard	34-35.5 cm
optional	37-38.5 cm
optional	43.5-45 cm

Castors	125 mm double-band 150 kg/castor
Steering castor	
standard	at foot end
optional	at head end
Lifting pole holder	
standard	outside one holder
optional	inside double holder
Pendant control	Autocontour
Caregiver pendant control	with total lock-outs
Design/Colours	see Völker Long Term Care beds brochure
Additional options	
Linen holder	

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN139315-03, 20 September 2016

Enhancing outcomes for patients and their caregivers:

Hill-Rom