

# Transport chair - Anatome

Makes transfer easier, enhances comfort and safety



## Features and benefits

- Adjusts to individual patients' shape, meaning improved comfort during transfers
- Intuitive design enhances patient comfort and is easy to operate for the caregiver
- Removable back and seat cushions allows for easier cleaning of the chair
- Height-adjustable armrests and retractable headrest enhance comfort and safety for patients
- Retractable footrest allows for a more comfortable seated position
- Multiple IV pole sockets and accessory hooks offer flexibility when placing accessories

## Transport chair - Anatome



Head support



Comfort seating pad



Overchair table



IV pole

### Technical Specifications

#### model SM647A

Height	1,345 mm
Seat height	510 mm
Width (with armrest)	680 mm
Depth (with footrest)	1,090 mm
Weight	45 kg
Maximum load	135 kg
Material	
Frame	Epoxy coated welded metal
Armrests	Polyurethane
Fabric	Coated polycotton blend
Castors	4 x 125 mm single-band with central brake

#### Fabrics

see coated fabric colours sheet

#### Accessories

- head support
- comfort seating pad
- overchair table
- telescopic IV pole (4 hooks)

#### model SM647B

Height	1,500 mm
Seat height (hydraulic)	495 - 715 mm
Width (with armrest)	680 mm
Depth (with footrest)	1,090 mm
Weight	58 kg
Maximum load	135 kg
Material	
Frame	Epoxy coated welded metal
Armrests	Polyurethane
Fabric	Coated polycotton blend
Castors	
standard	4 x 125 mm single-band with central brake
optional	4 x 125 mm double-band with central brake

#### Fabrics

see coated fabric colours sheet

#### Accessories

- head support
- comfort seating pad
- overchair table
- telescopic IV pole (4 hooks)

Hill-Rom is a leading global medical technology company with more than 10,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries.  
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN150323-01, 10 September 2015

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**