



## Roller-Free Roller Board: Better for Patients, Easier on Nurses and Inexpensive

The Hill-Rom Patient Transfer Board is a lightweight device designed to make lateral patient transfer easier for OR staff, and more comfortable for the patient. Its unique design is more comfortable for patients than traditional metal roller boards because of its soft foam core. It is also lightweight, for easy transportability and storage. In addition, Hill-Rom patient transfer boards are available in three sizes to accommodate different patient transfer needs throughout the hospital.

Ideal for use with:

- Transferring patients from OR tables to stretchers
- Transferring patients from stretchers to hospital beds

### KEY FEATURES AND BENEFITS

- **Comfortable** – Soft core with no hard metal rollers
- **Lightweight** – No lifting required, even with your heavy patients<sup>1</sup>
- **Simple** – Uses same or similar transfer technique as those used with metal roller boards
- **Fowler Position Transfer** – Long board allows for transfer in seated position
- **Foldable** – Long board folds for easy storage and set-up
- **Strap** – Carry this lightweight device on your shoulder and keep your hands free
- **Resilient and Sturdy** – Durable material protects against end cap breakage

<sup>1</sup> Max patient weight is 500 lbs

## About the Product

The new Hill-Rom Patient Transfer Boards use smooth, low friction surfaces to eliminate the need for hard metal rollers. The transfer board is lightweight, for easy handling and patient comfort. The device enables patient transfer in the seated or Fowler position. Both boards have a built-in strap for easy storage and transport. The Hill-Rom Patient Transfer Boards provide all of these advantages without requiring any change to transferring techniques used with traditional metal roller boards.

### Technical Specifications

#### Short – Narrow Board (P8300A06)

Weight	5.0 lbs (2.3 kg)
Dimensions	40" x 16.5" (100 cm x 42 cm)

#### Short – Wide Board (P8300A07)

Weight	4.4 lbs (2.0 kg)
Dimensions	32" x 20.5" (80 cm x 52 cm)

#### Long – Wide Board (P8300A08)

Weight	7.0 lbs (3.2 kg)
Dimensions	71.5" x 20.5" (180 cm x 52 cm)



SHORT-WIDE BOARD



SHORT-NARROW BOARD

LONG-WIDE BOARD

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

ORDER NUMBER 164450 rev 4 28-JAN-2014 ENG – US

Zinc Number US-AC-1113-0588

[www.hill-rom.com](http://www.hill-rom.com)

USA 800-445-3730  
CANADA 800-267-2337

Enhancing outcomes for  
patients and their caregivers:

