



Simply Easy to Use

The easy-to-use Liko® M220/M230 Mobile Lift is designed to meet lifting needs in home care environments with the durability and quality you've come to trust from Hill-Rom.

The Liko M220 Mobile Lift features an intuitive design intended to make the lift extraordinarily easy-to-use, even for caregivers new to lift operation. Coupled with a Liko sling, the lift creates a secure lifting experience for both patients and caregivers.

The Liko M230 lift features an electric lift mechanism with a powered base-width adjustment while the M220 Lift features a manual base width adjustment. The long lift arm on the Liko M220/M230 lift provides full access to the middle of the bed and its attractive, compact design works well in the home.



At a glance

- Electric high/low lifting motion
- Two-point sling bar makes sling attachment easy
- Gentle actuator speed for patient safety
- Quick Reference Guide reminds the caregiver how to operate the lift
- Low-profile 4.1" base leg clearance to fit easily under beds
- Smooth, non-marring casters are easy to clean and maneuver
- Manual base-width adjustment (M220)
- Powered base-width adjustment (M230)
- Long-lasting battery can be charged while on the lift
- Two-button (M220) and four-button (M230) hand control is easy-to-use
- Electrical and manual emergency lowering
- Emergency stop button quickly and safely stops all electrical functions



Electric high/low lifting motion

- Designed to provide safety and comfort for your patient and their caregiver

Dual-mode emergency lowering

- Electrical emergency lowering uses available battery power to safely lower resident
- Manual emergency lowering features knob with arrows showing direction to turn

Easy-to-use hand control

- Cable length allows caregiver to maintain contact with their patient while using the lift
- M230 four-button hand control for lifting and base width adjustment
- M220 two-button hand control for lifting

Quick reference guide

- Reminds caregivers how to operate lift

Long-lasting battery

- Charged while on the lift or easily removable for optional wall charging

Emergency stop button

- Quickly and safely stops electrical functions

Low battery indication

- Audible battery indication while pressing hand control button signals caregiver when battery is low

Long Lift Arm

- Provides access to center of bed
- Facilitates daily transfers of patients and can lift from the floor

Easy-to-use, standard two-point sling bar

- Lifts up to 400 lbs
- Compatible with Liko slings, scales and accessories

Gentle actuator speed for patient comfort and safety

Low-profile 4.1" leg clearance to fit easily under beds

- Front 3" twin wheels and rear 3" twin wheels with brakes are non-marring

Liko M230 electric base

- Allows caregiver to access their patient in a variety of lifting situations with the touch of a button

Liko M220 features a manual base with easy-to-use foot pedal for adjustments

Technical Specifications

Liko M220 (manual base)

Lifting capacity	400 lbs
Lifting interval	22.4" to 67.7"
Base-width adjustment	26.9" to 40"
Base length	45.7"
Wheels	Front 3" twin wheels, Rear 3" twin wheels with brakes
Material	Steel
Emergency lowering	Manual and Electrical
Battery	In series: 2x12 V 2.9 Ah valve-regulated lead-acid gel batteries
Battery charger	Built-in, 100-240 V AC, 50-60 Hz, max 400 mA
Total unit weight	90 lbs
Heaviest removable part	46 lbs

Technical Specifications

Liko M230 (electric base)

Lifting capacity	400 lbs
Lifting interval	22.4" to 67.7"
Base-width adjustment	26.9" to 40"
Base length	45.7"
Wheels	Front 3" twin wheels, Rear 3" twin wheels with brakes
Material	Steel
Emergency lowering	Manual and Electrical
Battery	In series: 2x12 V 2.9 Ah valve-regulated lead-acid gel batteries
Battery charger	Built-in, 100-240 V AC, 50-60 Hz, max 400 mA
Total unit weight	90 lbs
Heaviest removable part	46 lbs

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2012 Liko R&D AB. ALL RIGHTS RESERVED.

ORDER NUMBER 169228 rev 2 12/10/12 ENG - US

Zinc Number US-PA-1012-0028

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.®

USA 888-545-6671
Canada 800-267-2337

www.hill-rom.com