



Easy-to-move for everyday lifting.

The Viking® M and L lifts provide full-featured, mobile, lifting assistance that makes them ideal for most common lifting situations in many care environments. These durable aluminum lifts handle loads of up to 550 lbs (Viking® L) despite their low weight.

Key features at a glance

- Designed in aluminum – strong, durable, lightweight and environmentally friendly
- Suitable for most common lifting needs in all types of care environments
- Intelligent diagnostic system informs caregivers that the lift is safe to use and when service is required
- Optional armrests help support the patient during ambulation and make it easier for staff to maneuver the lift
- Light-weight, environmentally-friendly NiMH-battery that is easy to charge and easy to handle
- Can be equipped with quick-release hook system to quickly and easily change the sling bar or install a scale
- Dual control options to operate the lift – by hand control or control panel


The right lift for many situations.

The excellent weight-to-capacity ratio makes the Viking® lifts easy to maneuver, while features such as two lifting speeds and immediate charge status notification offer added convenience for caregivers.

Together with Hill-Rom's wide range of functional Liko® slings and accessories, the Viking M and Viking L lifts are ideal all-

around supports for those with high demands and wide-ranging needs.

The Viking patient lift is equipped with an intelligent diagnostic system that helps caregivers know when the lift is ready to use and when it is time for service. The lightweight, environmentally-friendly NiMH-battery is easy to charge and easy to handle.



	<u>Viking® M Lift</u>	<u>Viking® L Lift</u>
Product number	2040035	2040034
Maximum load	450 lbs	550 lbs
Weight	66.6 lbs	82.5 lbs
Lifting speed	1.7"/s	0.8"/s
Alternative lifting speed	1.5"/s	0.5"/s
Maximum lifting height	70.5"	70.5"
Base height	4.2"	6.5"
Width	27.2"	27.4"
Length	48.4"	53.5"
Material	Aluminum	Aluminum
Emergency lowering	Mechanical and electrical	Mechanical and electrical
Base-width adjustment	Electrical	Electrical

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2013 Liko R&D AB. ALL RIGHTS RESERVED.

ORDER NUMBER 179883 rev 1 E 05/09/2013 ENG – US

Zinc Number US-AC-1212-0079a

www.hill-rom.com

USA **800-445-3730**
CANADA **800-267-2337**