

Hill-Rom® OB-GYN Stretcher



Go Beyond Mobility™

With its exclusive Stow and Go® foot section, the Hill-Rom OB/GYN Stretcher is multi-functional for general use to pelvic exams. With a 700 lb weight capacity, integrated catch basin, and infinite labor positioning, our stretcher can be your stand-in birthing bed when you need it.

A rush to the operating room is not just about running down hallways and turning corners. It's about the minute details that contribute to a safe and comfortable experience for both you and your patient.

Our stretchers **provide more than mobility** for your patients:

- A safe and comfortable patient experience
- Simple and safe maneuvering
- Durability to withstand your demands

Key benefits at a glance:

A Safe and Comfortable Patient Experience

- Provide infinite positioning through foot and calf supports.
- Help provide comfort through Comfortline® foam mattress.
- Protect your patient's hands with bumper protected siderails.
- Provide a gentle patient ride with spring-loaded Steering Plus™ system that absorbs shock of driving over bumps and thresholds.

Simple and Safe Maneuvering

- Stow and Go® Foot Section requires no lifting.
- 360° maneuvering with Steering Plus™ system
- Zero Transfer Gap with OneStep® Tuck-Away Siderails

Durability to Withstand Your Demands

- 700 lb weight capacity
- Compatibility with most lifts and overbed tables
- Easy to clean smooth deck and siderail design
- Power washable

Go Beyond Mobility™

A safe and comfortable experience for both you and your patient starts on a Hill-Rom® stretcher.



Infinite Positioning Foot Supports



Integrated Catch Basin



Telescoping Calf Supports



Steering Plus™ System



Ergonomic Push Handles with exclusive IV Pole Holder



Permanent 3-Stage IV Pole

A SAFE & COMFORTABLE PATIENT EXPERIENCE	STANDARD FEATURES	OPTIONAL FEATURES
30" width	X	
Infinite Positioning Foot Supports	X	
Telescoping Calf Supports		X
Bumper Protected Siderails	X	
Trendelenburg/Reverse Trendelenburg	X	
Six Transport Strap Locations	X	
PRESSURE REDISTRIBUTION MATTRESS OPTIONS		
3" Comfortline® foam mattress	X	
SIMPLE & SAFE MANEUVERING		
Steering Plus™ System	X	
Integrated Catch Basin		X
OneStep® Tuck-Away Siderails	X	
Brake/Steer – head and foot	X	
Brake/Steer – all four sides		X
Dual-Sided hi/lo Hydraulic Control Pedals	X	
Ergonomic Push Handles	X	
Ergonomic Push Handles with Integrated IV Pole Transport Holder		X
3/4" Standard IV Pole Sockets	X	
Removable IV Pole – 2 Stage		X
Permanent IV Pole – 3 Stage (40 lb weight capacity)		X
Radiolucent Deck	X	
DURABILITY TO WITHSTAND YOUR DEMANDS		
700 lb weight capacity	X	
Full Perimeter Bumpering	X	
Smooth Deck Design	X	
Decals Options	X	
Bumper Color Options – choice of neutral, blue, red, teal or purple	X	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
197954 rev 1 29-JUN-2016 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Technical Specifications

Maximum weight limit	700 lbs (317.5 kg)
Low position (floor to top of sleep deck)	24.25" (61.6 cm)
High position (floor to top of sleep deck)	37" (94 cm)
Overall length	83" (210.8 cm)
Overall width (siderails up)	35.5" (90.2 cm)
Overall width (siderails stored)	32.25" (81.9 cm)
Siderail length	37" (94 cm)
Siderail height above sleep deck	13" (33 cm)
Comfortline® mattress size	30" x 72" (76.2 cm x 182.9 cm)
Maximum head elevation	90°
Maximum Trend/Reverse Trend	18°
Floor to base clearance	3.5" (8.9 cm)
Foot support inclination	70°
Caster size	8" (20.3 cm)

Warranty

Frames & welds for life of product

Three-year parts

Two-year mattress

One-year labor

Accessory Options — Please confirm accessory availability for model of stretcher ordered

Accessory IV pole	Liquid oxygen tank holder
Ankle stirrups	Oxygen tank holder
Armboard	Paper roll dispenser
Chart holder	Patient tray
Convertible footboard	Restraint straps
Detachable footboard	Siderail covers
Foot extender	Stretcher pillow
Infusion support system	Transport straps
IV transporter	Utility tray
Lateral X-ray cassette holder	

Enhancing outcomes for patients and their caregivers:

Hill-Rom