


Völker NP50 Monodensity

High-resilience foam mattress

Unverwechselbar
VÖLKER™



Efficiently conforming to your daily needs

- Durable and cost-effective patient support
- Turnable mattress
- High-resilience, monodensity foam
- Easy to clean and hygienic
- Adapts to most bed frames and conforms to all bed positions

Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Völker NP50 Monodensity

High-resilience foam mattress

Made of high-resilience foam, the Völker NP50 Monodensity mattress provides a high level of comfort for your patients. The Völker NP50 Monodensity mattress is adaptable to most bed frames and conforms to all bed positions, and the ability to turn the mattress promotes longevity.

The fully enclosed mattress has a waterproof, yet moisture vapour permeable, bacteriostatic and fungostatic polyurethane cover that is easy to clean and disinfect enhancing hygiene for patients and caregivers.



The Völker NP50 Monodensity offers 2 extended lengths

Standard size: 200 cm

+10 cm +20 cm

1

2

Technical Specifications

Mattress (customised sizes available):

Length	200 cm
Width	88 cm
Height	12 cm

Warranty:

Hill-Rom offers a 2 year warranty on the foam and the cover.
For details of the warranty provided please contact your local Hill-Rom representative.

Safe Working Load	250 kg
-------------------	--------

Recommended therapeutic weight	30 to 150 kg
--------------------------------	--------------

Disinfection:

Disinfect cover by wiping down with 1000 ppm chlorine-based solution

Conformity to standards:

EN ISO 10993-1:2009, EN ISO 10993-5:2009, EN ISO 10993-10:2009,
EN 597-1:1995, EN 597-2:1995, GPEM/CP - Guidelines D1-90 of February 1991,
SS 876 00 01 (NF FIRE 037)

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries.
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

The Völker NP 50 surface helps to prevent pressure ulcers in adult patients with no risk or at low risk, weighing between 30 kg and 150 kg with a minimum height of 1.46 m, in a range of healthcare environments.

Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 ZI du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN116308-02, 20 April 2017

Enhancing outcomes for
patients and their caregivers:

Hill-Rom