



# Hill-Rom® Enterprise Reporting

Protecting Patients by Anticipating Care<sup>SM</sup>

## Data Transparency that Drives Enhanced Patient Safety and Patient Satisfaction

Without reporting, medical device and hospital system data is limited in the impact it can provide.

Hospitals are filled with medical devices and systems that continuously generate data while supporting clinicians as they provide care to their patients. Reporting solutions play a crucial role in collecting all of that data and generating value.

Hill-Rom's Enterprise Reporting solution collects the information that matters and delivers it in a way that is relevant and actionable for caregivers and hospital leadership — all with the aim of improving care for patients, driving operational efficiencies, and enabling accountability.

**5,000** Beds, across multiple facilities, in one reporting database

**1,100+** Data points across nurse call, bed and real-time locating solutions

**20+** Built-in reports and dashboards

**Limitless** Options by connecting the database to third-party tools or using Hill-Rom Custom Report Services

# Hill-Rom® Enterprise Reporting

Protecting Patients by Anticipating Care<sup>SM</sup>

Hill-Rom's Enterprise Reporting solution can help clinical leadership assess a baseline of care, set goals for improvement, and impact the provision of care at the bedside.

## Server Requirements

Under 300 beds: Single Server

### Server 1: Report Model & Data Warehouse

Processor	1 Eight-Core 2.66 GHz
Memory	64 GB
OS	Windows Server 2012 R2
SQL	SQL Server 2016 Standard

Over 300 beds: Two Servers

### Server 1: Report Model

Processor	1 Eight-Core 2.66 GHz
Memory	64 GB RAM (up to 500 beds) Each additional 100 beds requires 5 GB RAM
OS	Windows Server 2012 R2
SQL	SQL Server 2016 Standard

### Server 2: Data Warehouse

Processor	1 Eight-Core 2.66 GHz
Memory	64 GB
OS	Windows Server 2012 R2
SQL	One of the following: <ul style="list-style-type: none"><li>• SQL Server 2016 Standard</li><li>• SQL Server 2016 Standard + CAL</li></ul>

## Opportunities to Impact Patient Safety and Patient Satisfaction:

- Improve compliance to rounding and falls protocols.
- Promote faster response times and more time with patients.
- Identify opportunities to proactively meet patient needs.
- Respond to patient complaints.

## Maximum Flexibility

When the need for a custom report arises, Hill-Rom's Enterprise Reporting database can be connected to a third-party tool of choice, like Tableau. This allows an informatics team to develop and deliver any custom reports that are needed.

Alternatively, Hill-Rom® Custom Report Services can provide a dedicated team with expertise and experience in working with nurse call, locating and smart bed data, to manage custom report creation from beginning to end.


Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
203604 rev 1 27-JUN-2017 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730  
Canada 800-267-2337

[www.hill-rom.com](http://www.hill-rom.com)



Enhancing outcomes for patients and their caregivers:

**Hill-Rom**